


CY5.5 DBCO

SKU: CCT-1046



Description


Cy5.5 DBCO is a bright and photostable near-IR probe that spectrally similar to Alexa Fluor® 680, DyLight® 680, and IRDye® 680 dye. The Cy5.5 DBCO is water-soluble, hydrophilic dye often a reagent of choice for assay where minimal non-specific binding and exceptional brightness is required. The fluorescence of Cy5.5 DBCO is pH insensitive from pH 4 to pH 10 and produces minimal autofluorescence of biological specimens in this region of the spectrum. Fluorescence of this long-wavelength Cyanine dye is not visible to the human eye but is readily detected by most imaging systems.

Cy5.5 DBCO reacts with azides via a copper-free “click chemistry” reaction to form a stable triazole and does not require Cu-catalyst or elevated temperatures. In application where the presence of copper is a concern Cy5.5 DBCO is an ideal alternative to copper requiring fluorescent alkynes.

Cy5.5 DBCO reagent is not suitable for staining intracellular components of fixed and permeabilized cells due to high backgrounds.

Abs/Em Spectra

For research use only. Not intended for animal or human therapeutic or diagnostic use.



Specifications

Unit Size	1 mg, 5 mg, 25 mg, 100 mg
Abs/Em Maxima	678/694 nm
Extinction Coefficient	190,000
Spectrally Similar Dyes	Alexa Fluor® 680, IRDye® 680RD, DyLight® 680
Molecular weight	1161.34
CAS	1857352-95-4
Solubility	Water, DMSO, DMF
Purity	>95% (HPLC)
Appearance	Blue solid
Storage Conditions	-20°C. Desiccate
Shipping Conditions	Ambient temperature

For research use only. Not intended for animal or human therapeutic or diagnostic use.