




DBCO ACID

SKU: CCT-1117



DESCRIPTION

DBCO Acid is an amine-reactive non-activated building block used for derivatizing amine-containing molecule in the presence of activators (e.g. EDC, or HATU) through a stable amide bond. Short spacer arm adds minimal mass to modified molecules (305.11 daltons).

DBCO reagents are the most commonly used substrates for copper-free click chemistry. DBCO compounds react with azide functionalized compounds or biomolecules without the need for a Cu(I) catalyst to result in a stable triazole linkage.

SPECIFICATIONS

CAS Number	1353016-70-2
Molecular Weight	305.11
Appearance	Slightly grey crystalline
Chemical Formula	C ₁₉ H ₁₅ N ₃ O ₃
Purity	>95% (HPLC)

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.




Unit Size	100 mg, 1000 mg, 5 g
Solubility	DMSO, DMF
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.