




DBCO AGAROSE

SKU: CCT-1034



DESCRIPTION

DBCO agarose is a 6% crosslinked agarose resin that is activated with DBCO functional groups for covalent capturing azide-tagged biomolecules utilizing Cu-free click reaction. This product is adequate to work in batch or column purification (low pressure).

SPECIFICATIONS

Activation Level	≥ 10 µmol DBCO groups per mL resin
Appearance	Off-white slurry
Bead Size	Spherical, 50-150 µm
Enrichment Target	Azide-modified proteins
Format	Batch or column purifications
Formulation	50% slurry
Preservative	20% ethanol
Unit Size	2 ml, 25 ml
Storage Instructions	2-8C
Support	6% Cross-linked agarose

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)


Shipping Conditions
Shipping Instructions

Ambient temperature
Ambient temperature

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[DBCO Agarose](#)

<https://vectorlabs.com/products/dbco-agarose/>