




## **DBC-O-AMINE**

**SKU:** QBD-11750



## **DESCRIPTION**

DBC-O-amine, product number QBD-11750, is an aza-dibenzocyclooctyne (DBC-O) product used for strain-promoted azide-alkyne cycloaddition (SPAAC, also known as “copper free click chemistry”). The DBC-O moiety reacts selectively with azides via SPAAC to form stable triazole linkages. The amine terminus of the molecule reacts with carboxylic acids and their active esters to form stable amide bonds. Insoluble in water, DBC-O-amine is slightly soluble to sparingly soluble in most organic solvents. Vector Laboratories recommends dimethylsulfoxide (DMSO), N,N'-dimethylformamide (DMF), or N,N'-dimethylacetamide (DMAC) as solvents for this product. DBC-O-amine is relatively insensitive to bases and can be used in the presence of aqueous bases. This product is acid sensitive. It should not be subjected to pH <5 during reaction or workup. The product should always be stored at -20°C when not in use.

The published discovery in 2004 of SPAAC, also known as “copper free click chemistry,” by Carolyn Bertozzi and coworkers improved the click chemistry toolbox. Toxic copper catalysts were eliminated, making the reaction truly bioorthogonal. By removing the copper(I) catalyst from the reaction, living cells can be labeled with reagents designed for SPAAC. Since the original publication, multiple improvements in strained cyclooctyne reagents resulted in the development of DBC-O as a nearly perfect azide partner for SPAAC.

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



To use DBCO-amine effectively in aqueous media or in cells, it should be conjugated to a water-soluble dPEG® product. Vector Laboratories offers numerous products for this purpose, including hydroxy acids (and their active esters), and protected amino acids. We also offer many DBCO-based dPEG® products already prepared for use.

## SPECIFICATIONS

<b>CAS Number</b>	N/A
<b>Molecular Weight</b>	276.33; single compound
<b>Chemical Formula</b>	C <sub>18</sub> H <sub>16</sub> N <sub>2</sub> O
<b>Purity</b>	> 97%
<b>Unit Size</b>	100 mg, 1000 mg
<b>Solubility</b>	Methylene Chloride, Methanol, Acetone, DMSO, or Acetonitrile.
<b>Spacers</b>	N/A
<b>Storage Instructions</b>	-20°C; Always let come to room temperature before opening; be careful to limit exposure to moisture and restore under an inert atmosphere; stock solutions can be prepared with dry solvent and kept for several days (freeze when not in use). dPEG® pegylation compounds are generally hygroscopic and should be treated as such. This will be less noticeable with liquids, but the solids will become tacky and difficult to manipulate, if care is not taken to minimize air exposure.
<b>Shipping Instructions</b>	Ambient

## DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

## GALLERY IMAGES


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



[www.vectorlabs.com](http://www.vectorlabs.com)

Email: [customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Telephone: [\(650\) 697-3600](tel:(650)697-3600)



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**

[DBC0-amine](#)

<https://vectorlabs.com/products/dbco-amine-2/>