




DBCO-dPEG®₁₂-MAL

SKU: QBD-10593



DESCRIPTION

DBCO-dPEG®₁₂-MAL, product number QBD-10593, is a medium-length (54 atoms), single molecular weight, discrete polyethylene glycol (dPEG®) crosslinking reagent. The dibenzylcyclooctyne (DBCO) moiety permits strain-promoted azide-alkyne cycloaddition (SPAAC), also known as “copper-free click chemistry.” The maleimidopropyl (MAL) functional group forms a thiol-ether linkage with sulfhydryl groups. The single molecular weight dPEG® spacer is hydrophilic and non-immunogenic. DBCO-dPEG®₁₂-MAL provides the benefits of PEG (increased water solubility, non-immunogenicity) without the analytical headaches of PEG dispersity.

The copper(I)-catalyzed azide-alkyne cycloaddition (CuAAC) discovered by K. Barry Sharpless and coworkers is one of the best-known click chemistry reactions, but its *in vivo* usability is limited by copper toxicity. Carolyn Bertozzi and coworkers developed SPAAC as a way to avoid copper's toxic effects on living cells using a highly strained cyclooctyne. Developments since the initial publication show that reagents based on DBCO provide excellent, highly selective reactivity with azides, forming a triazole product.

Maleimide groups are widely popular for their ability to react rapidly with sulfhydryl groups via a Michael addition to form thioether bonds.

Vector Laboratories' dPEG® products clearly differ from traditional polymeric polyethylene glycol (PEG). Traditional PEG is a dispersed polymer product with a Poisson distribution of chain lengths and molecular weights. The intractable mixture of chain lengths and molecular weights

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



in traditional PEG products creates significant analytical challenges. In contrast, Vector Laboratories' dPEG® products are single molecular weight compounds with discrete chain lengths. Using Vector Laboratories' products greatly simplifies the analysis of dPEG®-conjugated compounds.


SPECIFICATIONS

CAS Number	N/A
Molecular Weight	1027.161; single compound
Chemical Formula	C ₅₂ H ₇₄ N ₄ O ₁₇
Purity	> 97%
Unit Size	25mg, 100mg , 500mg
Solubility	Methylene Chloride, Methanol, DMF, DMSO, or DMAC.
Spacers	dPEG® Spacer is 54 atoms and 56.426 Å
Storage Instructions	-20°C; Always let come to room temperature before opening; be careful to limit exposure to moisture and restore under an inert atmosphere; stock solutions can be prepared with dry solvent and kept for several days (freeze when not in use). dPEG® pegylation compounds are generally hygroscopic and should be treated as such. This will be less noticeable with liquids, but the solids will become tacky and difficult to manipulate, if care is not taken to minimize air exposure.
Shipping Instructions	Ambient

DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[DBCO-dPEG®12-MAL](https://vectorlabs.com/products/dbco-dpeg12-mal/)

<https://vectorlabs.com/products/dbco-dpeg12-mal/>