




DBCO-MPEG, 30 KDA

SKU: CCT-A121



DESCRIPTION

DBCO-activated PEGylation reagents react with azides via a copper-free “click chemistry” reaction to form a stable triazole and does not require Cu-catalyst or elevated temperatures. In application where the presence of copper is a concern DBCO-activated PEGylation reagents is an ideal alternative to copper requiring fluorescent alkynes.

SPECIFICATIONS

CAS Number	N/A
Molecular Weight	30 kDa
Appearance	White crystalline
Unit Size	25mg, 100 mg, 1000 mg
Solubility	Water, DMSO, DMF, chloroform
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

DOCUMENTS


- [Safety Data Sheet](#)

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.