


DBCO-MPEG, 30 KDA

SKU: CCT-A121



Description

DBCO-activated PEGylation reagents react with azides via a copper-free “click chemistry” reaction to form a stable triazole and does not require Cu-catalyst or elevated temperatures. In application where the presence of copper is a concern DBCO-activated PEGylation reagents is an ideal alternative to copper requiring fluorescent alkynes.

Specifications

Unit Size	25mg, 100 mg, 1000 mg
Molecular weight	30 kDa
CAS	N/A
Solubility	Water, DMSO, DMF, chloroform
Appearance	White crystalline
Storage Conditions	-20°C. Desiccate
Shipping Conditions	Ambient temperature

For research use only. Not intended for animal or human therapeutic or diagnostic use.