




DDE AZIDE AGAROSE

SKU: CCT-1139



DESCRIPTION

Dde Azide Agarose is a 6% crosslinked agarose resin that is activated with azide groups through a cleavable Dde linker. Once biomolecules are covalently captured to the resin, beads can be washed with highest stringency virtually eliminating any non-specifically bound proteins. Treatment with 2% hydrazine yields a highly enriched population of proteins. This product is adequate to work in batch or column purification (low pressure).

SPECIFICATIONS

Activation Level	≥ 5 μmol Dde Azide groups per mL resin
Bead Size	Spherical, 50-150 μm
Enrichment Target	Alkyne-modified proteins
Format	Batch or column purifications
Formulation	50% slurry
Molecular Weight Left Behind	142.12
Unit Size	2 ml
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com


Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

Dde Azide Agarose

<https://vectorlabs.com/products/dde-azide-agarose/>