




FMOC-N-AMIDO-DPEG®₄-T-BOC-HYDRAZIDE

SKU: QBD-10043



DESCRIPTION

Fmoc-N-amido-dPEG®-t-boc-hydrazide, product number QBD-10043, is a heterobifunctional PEG modifier containing a short discrete PEG chain. It is useful in syntheses where an acid-labile hydrazone functionality is required.

SPECIFICATIONS

CAS Number	1263044-77-4
Molecular Weight	601.69; single compound
Chemical Formula	C ₃₁ H ₄₃ N ₃ O ₉
Purity	> 98%
Unit Size	100 mg, 1000 mg
Solubility	Methylene chloride, Acetonitrile, DMAC or DMSO.
Spacers	dPEG® Spacer is 18 atoms and 21.1 Å

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Storage Instructions


-20°C; Always let come to room temperature before opening; be careful to limit exposure to moisture and restore under an inert atmosphere; stock solutions can be prepared with dry solvent and kept for several days (freeze when not in use). dPEG® pegylation compounds are generally hygroscopic and should be treated as such. This will be less noticeable with liquids, but the solids will become tacky and difficult to manipulate, if care is not taken to minimize air exposure.

Shipping Instructions Ambient

DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.