



IMMPRINT PERMANENT MARKING PEN

SKU: H-6100



DESCRIPTION

The ImmPrint Permanent Marking Pen is a solvent-resistant permanent marking pen that is designed for writing on frosted glass microscope slides, tissue cassettes, and most hard surfaces.

Features:

- High density, fast-drying, black ink is resistant to most organic solvents encountered in histological applications such as alcohols and clearing agents
- Has a smooth writing tip that resists drying out

SPECIFICATIONS

Unit Size 5-Pen Set

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Applications Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.