




L-HOMOPROPARGYLGLYCINE (HPG)

SKU: CCT-1067



DESCRIPTION

L-Homopropargylglycine (HPG) is an amino acid analog of methionine that contains a very small modification, specifically an alkyne moiety that can be fed to cultured cells and incorporated into proteins during active protein synthesis. Detection utilizes the chemoselective ligation or “click ” reaction between an azide and an alkyne or cyclooctyne. The alkyne-modified protein is detected with either fluorescent azide or biotin azide. Detection sensitivity with these reagents in 1-D gels and Western blots is in the low femtomole range and compatible with downstream LC-MSMS and MALDI MS analysis.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



SPECIFICATIONS

CAS Number	942518-19-6
Molecular Weight	163.6
Molecular Formula	C ₆ H ₁₀ CINO ₂
Appearance	White crystalline
Chemical Formula	C ₆ H ₁₀ CINO ₂
Purity	>95% (H NMR)
Unit Size	25 mg, 100 mg, 1000 mg
Solubility	Water, DMSO, DMF
Storage Instructions	-20°C.
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

SELECTED REFERENCES


1. Couradeau, E., *et al.* (2019). Probing the active fraction of soil microbiomes using BONCAT-FACS *Nat Commun*, **10**, 2770-80. [[PubMed](#)]

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.