




MAL-DPEG® 4-(M-DPEG® 24)3

SKU: QBD-10456



DESCRIPTION

MAL-dPEG® 4-(m-dPEG® 24)3, product number QBD-10456, is a sulfhydryl-reactive, branched polyethylene glycol (PEG) product built around a tris core that is useful for modifying surfaces and biomolecules containing free thiols. Vector Laboratories, Inc. manufactures this monodispersed PEG product with a single molecular weight and discrete chain length and sells it under the dPEG® tradename. MAL-dPEG® 4-(m-dPEG® 24)3 increases the hydrophilicity of conjugated molecules; reduces non-specific, hydrophobic interactions; and decreases or eliminates the antigenicity of conjugates. Also, because it is a single molecule, the analysis of conjugates made with MAL-dPEG® 4-(m-dPEG® 24)3 is simplified compared to the analysis of conjugates prepared with traditional dispersed PEG products.

SPECIFICATIONS

CAS Number	1334178-05-0
Molecular Weight	3946.64; single compound
Chemical Formula	C ₁₇₈ H ₃₄₆ N ₆ O ₈₆
Purity	> 96%
Unit Size	100 mg, 1000 mg
Solubility	Methylene chloride, DMAC or DMSO.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Spacers

dPEG® Spacer is 103 atoms and 108.3 Å, Avg.

Storage Instructions


-20°C; Always let come to room temperature before opening; be careful to limit exposure to moisture and restore under an inert atmosphere; stock solutions can be prepared with dry solvent and kept for several days (freeze when not in use). dPEG® pegylation compounds are generally hygroscopic and should be treated as such. This will be less noticeable with liquids, but the solids will become tacky and difficult to manipulate, if care is not taken to minimize air exposure.

Shipping Instructions Ambient

DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.