




METHYLTETRAZINE-MPEG, 10 KDA

SKU: CCT-1091



DESCRIPTION

Methyltetrazine-activated PEGylation reagent that reacts with TCO-containing compounds via Inverse-Electron-Demand Diels-Alder reaction to form a stable covalent bond and does not require Cu-catalyst or elevated temperatures.

SPECIFICATIONS

CAS Number	N/A
Molecular Weight	10 kDa
Appearance	Pink solid
Unit Size	25 mg, 100 mg, 1000 mg
Solubility	Water, DMSO, DMF, DCM, Chloroform
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

DOCUMENTS


- [Safety Data Sheet](#)

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.