


METHYLTETRAZINE-MPEG, 5 KDA

SKU: CCT-1090



Description

Methyltetrazine-activated PEGylation reagent that reacts with TCO-containing compounds via Inverse-Electron-Demand Diels-Alder reaction to form a stable covalent bond and does not require Cu-catalyst or elevated temperatures.

Specifications

Unit Size	25 mg, 100 mg, 1000 mg
Molecular weight	5 kDa
CAS	N/A
Solubility	Water, DMSO, DMF, DCM, Chloroform
Appearance	Pink solid
Storage Conditions	-20°C. Desiccate
Shipping Conditions	Ambient temperature

For research use only. Not intended for animal or human therapeutic or diagnostic use.