




METHYLTETRAZINE-PEG4-ALKYNE

SKU: CCT-1013



DESCRIPTION

Methyltetrazine-PEG4-Alkyne is a heterobifunctional, azide-to-TCO crosslinker with terminal alkyne and methyltetrazine groups connected through a PEG4 spacer. Often used to convert a readily available azide groups into methyltetrazine moiety through CuAAC reaction.

SPECIFICATIONS

CAS Number	N/A
Molecular Weight	487.56
Appearance	Red amorphous solid to red oil
Chemical Formula	C ₂₄ H ₃₃ N ₅ O ₆
Purity	>95% (HPLC)
Unit Size	10 mg, 25 mg, 100 mg
Solubility	DMSO, DMF, DCM, THF, Chloroform
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.