



MODIFICATION BUFFER (10X)

SKU: S-4000-005



DESCRIPTION

Modification Buffer is an amine-free, high buffering capacity solution which is optimized for the reaction of NHS esters and TFP esters with amino groups. Proteins, oligonucleotides, and other biomolecules are buffer exchanged into Modification Buffer prior to modification with S-HyNic and S-4FB to ensure efficient coupling of these linkers.

SPECIFICATIONS

Unit Size

5 x 1.5 mL

Storage Instructions

2° - 8°C - Do Not Freeze

Applications

In Situ Proximity Ligation

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



LEGAL NOTICE

Products are for research use only, not for use in diagnostic or therapeutic procedures or for use in humans. Products are not for resale without express written permission of Seller. No license under any patent or other intellectual property right of Seller or its licensors is granted or implied by the purchase unless otherwise provided in writing.

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.