



MOUSE IGG, CONTROL ANTIBODY

SKU: I-2000-1



DESCRIPTION

This IgG preparation is intended for use as a control for primary antibodies made in mouse. This antibody has been purified from pooled serum of healthy adult animals and contains a spectrum of the IgG subclasses present in serum. The control antibody should be applied to the tissue section at the same concentration as the primary antibody to indicate whether staining is specific for the antigen or is nonspecific adsorption of primary antibody to tissue sites.

SPECIFICATIONS

Unit Size 1 mg

Storage Instructions 2-8 °C; Store frozen for long term storage

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Usage Summary

Reconstitute by adding 0.5 ml water. The resulting solution will have the following composition: 10 mM phosphate, 0.15 M NaCl, pH 7.8, 0.08% sodium azide, 20 mg/ml sucrose. ALLOW THE PROTEIN TO DISSOLVE WITHOUT VIGOROUS AGITATION.

Applications

Immunohistochemistry / Immunocytochemistry,
Immunofluorescence, Blotting Applications

Conjugate

Unconjugated

Reactive Species

Mouse

Host Species

Mouse

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.