




MPS-EDA.TFA

SKU: QBD-10177



DESCRIPTION

MPS-EDA•TFA, product number QBD-10177, is the trifluoroacetic acid salt of N-(2-aminoethyl)-3-(2,5-dioxo-2,5-dihydro-1H-pyrrol-1-yl)propanamide. This product crosslinks a sulfhydryl group and a carboxylic acid. The spacer length from the maleimide reactive site to the terminal nitrogen is ten (10) atoms long (10.7 Å). The product is a highly water-soluble white solid. It is extremely stable under recommended storage conditions. Although this product does not contain a dPEG® linker or spacer, it can be conjugated to the carboxylic acid termini of various dPEG® products to form maleimide-functionalized constructs.

SPECIFICATIONS

CAS Number	11550-02-00
Molecular Weight	325.24; single compound
Chemical Formula	$\text{C}_{11}\text{H}_{14}\text{F}_3\text{N}_3\text{O}_5$
Purity	> 98%
Unit Size	100 mg, 1000 mg
Solubility	DMAC, DMSO or water.
Spacers	Spacer is 10 atoms and 10.7 Å

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Storage Instructions


-20°C; Always let come to room temperature before opening; be careful to limit exposure to moisture and restore under an inert atmosphere; stock solutions can be prepared with dry solvent and kept for several days (freeze when not in use). dPEG® pegylation compounds are generally hygroscopic and should be treated as such. This will be less noticeable with liquids, but the solids will become tacky and difficult to manipulate, if care is not taken to minimize air exposure.

Shipping Instructions Ambient

DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.