



# **NORMAL CHICKEN SERUM BLOCKING SOLUTION**

**SKU:** S-3000-20



---

## **DESCRIPTION**

Normal Chicken Serum can be used for blocking non-specific binding or as an antibody diluent.

### **Features:**

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Supplied undiluted with 0.08% sodium azide as preservative

## **SPECIFICATIONS**

**Blocking Action** Non-Specific Protein Blocking

**Format** Concentrate

**Unit Size** 20 ml

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



**Applications** — Immunohistochemistry / Immunocytochemistry, Immunofluorescence

## CITATIONS



Powered by Bioz © 2023 See more details on Bioz

## DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**