



NORMAL CHICKEN SERUM BLOCKING SOLUTION

SKU: S-3000-20



DESCRIPTION

Normal Chicken Serum can be used for blocking non-specific binding or as an antibody diluent.

Features:

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Supplied undiluted with 0.08% sodium azide as preservative

SPECIFICATIONS

Blocking Action Non-Specific Protein Blocking

Format Concentrate

Unit Size 20 ml

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Applications — Immunohistochemistry / Immunocytochemistry, Immunofluorescence

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[Normal Chicken Serum Blocking Solution](https://vectorlabs.com/products/normal-chicken-serum-blocking-solution/)

<https://vectorlabs.com/products/normal-chicken-serum-blocking-solution/>