



NORMAL GOAT SERUM BLOCKING SOLUTION

SKU: S-1000-20



DESCRIPTION

Normal goat serum can be used for blocking non-specific binding or as an antibody diluent.

Features:

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Supplied undiluted with 0.08% sodium azide as preservative

SPECIFICATIONS

Blocking Action Non-Specific Protein Blocking

Format Concentrate

Unit Size 20 ml

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Applications Immunohistochemistry / Immunocytochemistry, Immunofluorescence

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.