



NORMAL HORSE SERUM BLOCKING SOLUTION

SKU: S-2000-20



DESCRIPTION

Normal horse serum can be used for blocking non-specific binding or as an antibody diluent.

Features:

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Supplied undiluted with 0.08% sodium azide as preservative

SPECIFICATIONS

| | |
|------------------------|-------------------------------|
| Blocking Action | Non-Specific Protein Blocking |
| Format | Concentrate |
| Unit Size | 20 ml |

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Applications — Immunohistochemistry / Immunocytochemistry, Immunofluorescence

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.