



# **NORMAL RABBIT SERUM BLOCKING SOLUTION**

**SKU:** S-5000-20



---

## **DESCRIPTION**

Normal rabbit serum can be used for blocking non-specific binding or as an antibody diluent.

### **Features:**

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Supplied undiluted with 0.08% sodium azide as preservative

## **SPECIFICATIONS**

<b>Blocking Action</b>	Non-Specific Protein Blocking
<b>Format</b>	Concentrate
<b>Unit Size</b>	20 ml

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



**Applications** — Immunohistochemistry / Immunocytochemistry, Immunofluorescence

## CITATIONS



Powered by Bioz © 2023 See more details on Bioz

## DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**