




PC ALKYNE AGAROSE

SKU: CCT-1142



DESCRIPTION

PC Alkyne Agarose is a 6% crosslinked agarose resin that is activated with alkyne groups through a photocleavable linker. Once biomolecules are covalently captures to the resin, beads can be washed with highest stringency virtually eliminating any non-specifically bound proteins. Irradiation with a hand-held UV-Vis lamp allows for reagent-less release of captured molecules. This product is adequate to work in batch or column purification (low pressure).

SPECIFICATIONS

Activation Level	≥ 5 μmol PC Alkyne groups per mL resin
Bead Size	Spherical, 50-150 μm
Enrichment Target	Azide-modified proteins
Format	Batch or column purifications
Formulation	50% slurry
Molecular Weight Left Behind	55.04
Unit Size	2 ml
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.