




Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:6506973600)

SKU: CCT-1162



Photocleavable Alkyne-NHS ester reagent contains an amine reactive group linked to terminal alkyne moiety through a spacer arm containing a photocleavable moiety. Molecules joined through photocleavable linker can be efficiently photoreleased, typically >90% in 5-25 minutes using an inexpensive, near-UV, low intensity lamp (e.g. 365 nm lamp at 1-5 mW/cm²).

CAS Number	N/A
Molecular Weight	656.63
Appearance	Yellow oil
Chemical Formula	C29H29N3O14
Purity	>90% (HPLC)
Unit Size	10 mg, 25 mg, 100 mg
Solubility	DMSO, DMF, DCM, THF, Chloroform
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Frozen
Shipping Instructions	Frozen

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

PC Alkyne-NHS Ester


<https://vectorlabs.com/products/pc-alkyne-nhs-ester/>



DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.