




PC DBCO BIOTIN

SKU: CCT-1120



DESCRIPTION

PC Biotin DBCO is an azide-reactive photocleavable biotin probe that allows for reagent-free release of the captured biomolecules from streptavidin. This reagent contains a biotin moiety linked to DBCO moiety through a spacer arm containing a photocleavable moiety. Captured biomolecules can be efficiently photoreleased, typically >90% in 5-25 minutes using an inexpensive, near-UV, low intensity lamp (e.g. 365 nm lamp at 1-5 mW/cm²).

SPECIFICATIONS

CAS Number	N/A
Molecular Weight	1002.14
Appearance	Yellow amorphous solid
Chemical Formula	C ₅₀ H ₆₁ N ₇ O ₁₃ S
Purity	>95% (HPLC)
Unit Size	1 mg, 5 mg, 25 mg
Solubility	DMSO, DMF, DCM, THF, Chloroform
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Frozen
Shipping Instructions	Frozen

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com


Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:6506973600)

DOCUMENTS

- [Safety Data Sheet](#)
- [Certificate of Analysis](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

PC DBCO Biotin

<https://vectorlabs.com/products/pc-dbc0-biotin/>