



# GOAT ANTI-RABBIT IGG ANTIBODY (H+L), PEROXIDASE

**SKU:** PI-1000-1



---

## DESCRIPTION

The HRP Goat Anti-Rabbit IgG Antibody binds to rabbit primary antibodies for tissue staining, ELISAs, and blots. It can be detected using a horseradish peroxidase substrate.

### Features:

- Recognizes both heavy and light chains (H+L)
- Cross-linked with the highest specific activity horseradish peroxidase
- Supplied in solution

## SPECIFICATIONS

<b>Format</b>	Concentrate
<b>Formulation</b>	10 mM HEPES, 0.15 M NaCl, pH 7.5, preservative
<b>Unit Size</b>	1 mg

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



**Storage Instructions** 2-8 °C, DO NOT FREEZE

**Usage Summary**

Recommended diluent for microtitration plate assays: 10 mM phosphate, 0.15 M NaCl, pH 7.5, 0.1% Tween 20. Working concentration of HRP-anti rabbit IgG is 0.5-5.0 µg/ml. Do not use sodium azide with this product since it is an inhibitor of horseradish peroxidase.

**Applications**

Immunohistochemistry / Immunocytochemistry, In situ hybridization, Blotting Applications, Elispot, ELISAs

**Target Species**

Rabbit

**Concentration**

1 mg active conjugate/ml

**Conjugate**

Peroxidase

**Safety Title**

WARNING: Cancer -[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Reactive Species**

Goat

**Source Species**

Rabbit

**Host Species**

Goat

## TECHNICAL INFORMATION

Our high affinity Goat Anti-Rabbit IgG Antibody is purified by affinity chromatography and adsorbed to remove any cross-reactivities that are likely to interfere with specific labeling. A proprietary procedure for conjugation allows the maximum preservation of enzyme activity and antibody specificity.

## CITATIONS



Powered by Bioz © 2023 See more details on Bioz

## DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



[www.vectorlabs.com](http://www.vectorlabs.com)

Email: [customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**

[Goat Anti-Rabbit IgG Antibody \(H+L\), Peroxidase](https://vectorlabs.com/products/peroxidase-goat-anti-rabbit-igg/)

<https://vectorlabs.com/products/peroxidase-goat-anti-rabbit-igg/>