




SPECIFICATIONS

CAS Number	N/A
Molecular Weight	860.90; single compound
Chemical Formula	C ₃₉ H ₆₀ N ₂ O ₁₉
Purity	> 97%
Unit Size	100mg, 1000mg
Solubility	Methylene chloride or Acetonitrile.
Spacers	dPEG® Spacer is 41 atoms and 46.3 Å
Storage Instructions	-20°C; Always let come to room temperature before opening; be careful to limit exposure to moisture and restore under an inert atmosphere; stock solutions can be prepared with dry solvent and kept for several days (freeze when not in use). dPEG® pegylation compounds are generally hygroscopic and should be treated as such. This will be less noticeable with liquids, but the solids will become tacky and difficult to manipulate, if care is not taken to minimize air exposure.
Shipping Instructions	Ambient

DOCUMENTS

- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.