




PICOLYL AZIDE AGAROSE

SKU: CCT-1408



DESCRIPTION

Picolyl Azide Agarose is a 6% crosslinked agarose resin that is activated with picolyl azide functional groups for covalent capturing alkyne-tagged biomolecules. Picolyl Azide incorporates a copper-chelating motif to raise the effective concentration of Cu(I) at the reaction site to boosted the efficiency of the CuAAC reaction resulting in a faster and more biocompatible CuAAC labeling. Up to 40-fold increase of signal intensity, compared to conventional azides was reported (see Selected References). This product is adequate to work in batch or column purification (low pressure).

SPECIFICATIONS

Activation Level	≥ 8 μmol Picolyl Azide groups per mL resin
Appearance	Off-white slurry
Bead Size	Spherical, 50-150 μm
Enrichment Target	Alkyne-tagged proteins
Format	Batch or column purifications
Formulation	50% slurry
Preservative	20% ethanol
Unit Size	2 ml, 25 ml
Storage Instructions	2-8C
Support	6% Cross-linked agarose
Shipping Conditions	Ambient temperature

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.




Shipping Instructions

Ambient temperature

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.