



Concanavalin A (Con A), Unconjugated

L-1000-500

[Product Images](#)



Short Description

Con A recognizes α -linked mannose present as part of a “core oligosaccharide” in many serum and membrane glycoproteins.

Additional Information

Unit Size	500 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, Blotting Applications, Glycobiology, Mitogenic Stimulation
Recommended Usage	Although many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM CaCl ₂ is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used.
Conjugate	Unconjugated
Sugar Specificity	Mannose, Glucose

