



---

# Normal Horse Serum Blocking Solution

## S-2000-20

[Product Images](#)

---



## Short Description

---

Normal horse serum can be used for blocking non-specific binding or as an antibody diluent.

### Features:

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Supplied undiluted with 0.08% sodium azide as preservative

## Additional Information

---

Unit Size	20 ml
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence
Blocking Action	Non-Specific Protein Blocking
Format	Concentrate

