



---

# N-acetylgalactosamine

## S-9001

[Product Images](#)

---



## Short Description

---

### Features:

- Provided for use as an inhibitor of lectin-conjugate binding or for eluting glycoproteins or other glucoconjugates from columns of agarose-lectins
- Supplied as dry powder
- If reconstituted in 5 ml of water or buffer, the resulting sugar solution is 100 mM

## Additional Information

---

Unit Size	111 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, Glycobiology, Affinity Chromatography
Recommended Usage	Reconstitute the sugar in 5 ml of distilled or deionized water or in 5 ml of a buffer (e.g. 10 mM HEPES buffered saline, pH 7.5-pH 8.0) or salt solution compatible with the lectin under study. If reconstituted in 5 ml, this sugar will be 100 mM.
Recommended Storage	Reconstituted sugar solutions should be stored frozen or stored at 4 °C using sodium azide (0.08%) as a preservative.

