



Anti-Streptavidin Antibody, Unconjugated

SP-4000-1

[Product Images](#)



Short Description

Anti-Streptavidin antibodies have been widely used as amplifying reagents in immunohistochemistry, *in situ* hybridization, microarray assays, ELISAs, blots, and many other applications. Our antibodies to streptavidin are produced in goats using our highly purified streptavidin and isolated by affinity chromatography. Anti-Streptavidin does not bind avidin and Anti-Avidin does not recognize streptavidin. These antibodies provide opportunities to significantly amplify signals in many applications.

Additional Information

Unit Size	1 mg
Applications	Immunofluorescence, In situ hybridization
Concentration	2 mg active IgG/ml
Recommended Storage	2-8 °C; Store frozen for long term storage
Solution	10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose.
Recommended Usage	The recommended concentration range for use is 2-5 µg/ml.
Detection Target	Streptavidin
Conjugate	Unconjugated
Host Species	Goat
Format	Concentrate

