



PROTEIN A, BIOTINYLATED (S. AUREUS, COWAN STRAIN)

SKU: B-2001-1



DESCRIPTION

Biotinylated Protein A (B-2001), isolated from *S. aureus*, is a 41 kDa protein which binds to the Fc portion of most immunoglobulins. This interaction is neither species specific nor isotype specific, as Protein A has been shown to bind IgG from a variety of species. Protein A can also interact with IgA, IgM and subclasses of IgG from different species.

SPECIFICATIONS

Format	Concentrate
Unit Size	1 mg
Storage Instructions	2-8 °C after reconstitution; -20 °C is recommended.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Usage Summary

Reconstitute by adding 0.5 ml water. The resulting solution will have the following composition: 0.15 M NaCl, 10 mM HEPES, pH 8.0 The recommended concentration range for use is 2-5 µg/ml. For long term storage, an immunohistochemical grade, globulin-free bovine serum albumin (0.1%, w/v) may be added to the biotinylated Protein A solution prior to freezing. For short-term usage, bovine serum albumin may also be added and samples stored at 4 °C in the presence of 0.05% sodium azide.

Applications

Immunohistochemistry / Immunocytochemistry,
Immunofluorescence, Blotting Applications, Elispot, ELISAs

Conjugate

Biotinylated

TECHNICAL INFORMATION

Biotinylated Protein A is useful to detect or localize immunoglobulins on cell surfaces or in tissue sections. Used with the appropriate avidin or streptavidin conjugates, immunoglobulins can be studied with brightfield, fluorescence or electron microscopy. Biotinylated Protein A can also be used to detect positive clones in monoclonal antibody production by an avidin or streptavidin enzyme conjugate.

Biotinylated Protein A is used in many systems as a substitute for a biotinylated second antibody. However, in some applications, Biotinylated Protein A may bind to endogenous immunoglobulin in the tissue producing background. In addition, Biotinylated Protein A may not bind to the species or class of primary antibody being employed. For routine immunohistochemical staining we recommend using biotinylated anti-immunoglobulins.

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[Protein A, Biotinylated \(S. Aureus, Cowan strain\)](#)

<https://vectorlabs.com/products/protein-a-biotinylated-s-aureus/>