

# Q-BRIGHT® ECL PLUS CHEMILUMINESCENT DETECTION KIT

**SKU:** QBD-60012



Q-Bright® ECL Plus chemiluminescent detection kit is a luminol-based substrate for horseradish peroxidase (HRP) providing a more robust signal over the original ECL reagents. Q-Bright® ECL Plus provides picogram to high-femtogram levels of detection in western blot, ELISA and other chemiluminescent-based applications.

## Specifications

<b>Unit Size</b>	20 ml, 100 ml, 200 ml
<b>Format</b>	Solution
<b>Storage</b>	Store at 4°C. Product may ship at ambient temperature.
<b>Country of Origin</b>	United States of America
<b>Sensitivity</b>	Picogram to high-femtogram

---

**For research use only. Not intended for animal or human therapeutic or diagnostic use.**