

# Q-BRIGHT® ENDURE CHEMILUMINESCENT DETECTION KIT

SKU: QBD-60013



Q-Bright® Endure chemiluminescent detection reagent is a luminol-based substrate for horseradish peroxidase (HRP), providing a light output for up to 24h. Q-Bright® Endure provides much greater sensitivity compared to the classical ECL reagents and is superior to many of the available commercial reagents providing long-lasting, high to mid-femtogram levels of detection in western blot, ELISA and other chemiluminescent-based applications.

## Specifications

<b>Unit Size</b>	20 ml, 100 ml, 200 ml
<b>Format</b>	Solution
<b>Storage</b>	Store at 4°C.
<b>Country of Origin</b>	United States of America
<b>Sensitivity</b>	high to mid-femtogram

---

**For research use only. Not intended for animal or human therapeutic or diagnostic use.**