

Q-BRIGHT® ENDURE PLUS CHEMILUMINESCENT DETECTION KIT

SKU: QBD-60014



Q-Bright® Endure Plus chemiluminescent detection reagent is a luminol-based substrate for horseradish peroxidase (HRP), providing a light output for up to 24h. Q-Bright® Endure Plus provides 2X the signal of Q-Bright® Endure without compromising on signal duration, providing long-lasting, mid-femtogram levels of detection in western blot, ELISA and other chemiluminescent-based applications.

Specifications

Unit Size	20 ml, 100 ml, 200 ml
Format	Solution
Storage	Store at 4°C.
Country of Origin	United States of America
Sensitivity	mid-femtogram

For research use only. Not intended for animal or human therapeutic or diagnostic use.