



RABBIT IGG, CONTROL ANTIBODY

SKU: I-1000-5



DESCRIPTION

This IgG preparation is intended for use as a control for primary antibodies made in rabbit. This antibody has been purified from pooled serum of healthy adult animals and contains a spectrum of the IgG subclasses present in serum.

SPECIFICATIONS

Unit Size	5 mg
Storage Instructions	2-8 °C; Store frozen for long term storage
Usage Summary	Reconstitute by adding 1 ml water. The resulting solution will have the following composition: 10 mM phosphate, 0.15 M NaCl, pH 7.8, 0.08% sodium azide, 20 mg/ml sucrose. ALLOW THE PROTEIN TO DISSOLVE WITHOUT VIGOROUS AGITATION.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, Blotting Applications
Conjugate	Unconjugated
Reactive Species	Rabbit
Host Species	Rabbit

TECHNICAL INFORMATION

The control antibody should be applied to the tissue section at the same concentration as the primary antibody to indicate whether staining is specific for the antigen or is nonspecific adsorption of primary antibody to tissue sites.

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.