




SY-21 DBCO

SKU: FP-1715



DESCRIPTION

SY-21 DBCO (QSY®-21 equivalent) is a nonfluorescent acceptor dye often used for preparation of peptide and oligonucleotide FRET probes. A broad and intense quenching range from 580 nm to 680 nm makes it useful as an acceptor in fluorescence resonance energy transfer (FRET) applications in conjunction with fluorescent dyes at 580 nm to 680 nm, such as Fluorescein, Alexa Fluor® 568, 594, 633, 647, TAMRA, ROX, Texas Red and Cy5.

SY-21 DBCO reacts with azides via a copper-free “click chemistry” reaction to form a stable triazole and does not require Cu-catalyst or elevated temperatures. In application where the presence of copper is a concern SY-21 DBCO is an ideal alternative to copper requiring fluorescent alkynes.


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



SPECIFICATIONS

Molecular Weight	1091.27
Extinction Coefficient	90,000 cm ⁻¹ M ⁻¹
Reactivity	Azide
Unit Size	5 mg, 25 mg, 100 mg
Solubility	DMSO, DMF
Storage Instructions	-20°C.
Spectrally Similar Dyes	Alexa Fluor® 568, 594, 633, 647, TAMRA, ROX, Texas Red, Cy5.
Excitation/Emission Maximum	661 nm
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

ABS/EM SPECTRA




DOCUMENTS

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



- [Safety Data Sheet](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.