




TAMRA BIOTIN ALKYNE

SKU: CCT-1366



DESCRIPTION

Whilst biotin is a highly effective affinity label, fluorescent labels provide a more sensitive, quantitative, and convenient method for visualizing proteins. Trifunctional click chemistry probes that incorporate a ligation handle, a biotin and a fluorophore straightforwardly extended to combine the complementary benefits of both types of label. Dual-label TAMRA Biotin Azide probe can be readily incorporated into alkyne-tagged biomolecules through CuAAC ligation. For example, dual-labeled newly synthesized proteins, glycans, lipids, or DNA/RNA can be readily isolated by affinity purification and visualized by in-gel fluorescence.

Another very important advantage of dual label probes over regular biotin probes is built-in control. Each step of enrichment process can easily followed either by UV-Vis (550 nm) or by more sensitive fluorescence spectroscopy. After elution from streptavidin beads target proteins containing TAMRA label can be easily distinguished from non-specifically bound proteins and endogenously biotinylated proteins.

This new dual label probe should be useful to expand the new applications of standard biotin-streptavidin technology.

SPECIFICATIONS

CAS Number

N/A

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.




Molecular Weight	1188.42
Appearance	Red solid
Chemical Formula	$C_{60}H_{83}N_8O_{15}S$
Unit Size	1 mg, 5 mg, 25 mg
Solubility	DMSO, DMF
Storage Instructions	-20°C.
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.