




TCO-PEG4-ACID

SKU: CCT-1088



DESCRIPTION

An amine-reactive, non-activated building block with enhanced solubility in aqueous media used for derivatizing amine-containing molecule in the presence of activators (e.g. EDC, or HATU) through a stable amide bond.

SPECIFICATIONS

CAS Number	N/A
Molecular Weight	417.49
Appearance	Colorless to slightly yellow oil
Chemical Formula	C ₂₀ H ₃₅ NO ₈
Purity	>95% (HPLC)
Unit Size	25 mg, 100 mg, 500 mg
Solubility	DMSO, DMF, DCM, THF, Chloroform
Storage Instructions	-20°C. Desiccate
Shipping Conditions	Dry ice
Shipping Instructions	Dry ice

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

