


THIOLINKER-DBCO

SKU: CCT-1096



Description

ThioLinker-DBCO is a water-soluble bis-alkylating labeling reagent that enables site-specific incorporation of DBCO moiety onto antibodies to form well-defined conjugates. This approach involves two steps: (1) disulfide reduction to release the two cysteine thiols and (2) re-forming the disulfide by bis-alkylation via a three-carbon bridge to which DBCO moiety is covalently attached. During this process, irreversible denaturation of the protein, disrupting either their tertiary structure or their biological activity does not occur.

The covalent, site-specific conjugation of PEG to a polyhistidine tag (His-tag) on a protein with bis-alkylating reagents was also reported in literature (see Selected References section).

Specifications

Unit Size	10 mg, 25 mg, 100 mg
Molecular weight	1140.39
Chemical composition	C ₅₉ H ₆₉ N ₃ O ₁₄ S ₃
CAS	N/A
Solubility	Water, DMSO, DMF DCM, THF, Chloroform
Purity	>95% (HPLC)
Appearance	White to grey amorphous solid
Storage Conditions	-20°C. Desiccate
Shipping Conditions	Ambient temperature

For research use only. Not intended for animal or human therapeutic or diagnostic use.

For research use only. Not intended for animal or human therapeutic or diagnostic use.