



ANTI-STREPTAVIDIN ANTIBODY, UNCONJUGATED

SKU: SP-4000-1



DESCRIPTION

Anti-Streptavidin antibodies have been widely used as amplifying reagents in immunohistochemistry, *in situ* hybridization, microarray assays, ELISAs, blots, and many other applications. Our antibodies to streptavidin are produced in goats using our highly purified streptavidin and isolated by affinity chromatography. Anti-Streptavidin does not bind avidin and Anti-Avidin does not recognize streptavidin. These antibodies provide opportunities to significantly amplify signals in many applications.

SPECIFICATIONS

Format	Concentrate
Formulation	10 mM HEPES, pH 7.5, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose.
Unit Size	1 mg

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Storage Instructions 2-8 °C; Store frozen for long term storage

Usage Summary The recommended concentration range for use is 2-5 µg/ml.

Applications Immunofluorescence, In situ hybridization

Concentration 2 mg active IgG/ml

Conjugate Unconjugated

Reactive Species Goat

Host Species Goat

TECHNICAL INFORMATION

These affinity purified antibodies are also available biotin or fluorescein-labeled.

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[Anti-Streptavidin Antibody, Unconjugated](https://vectorlabs.com/products/unconjugated-anti-streptavidin/)

<https://vectorlabs.com/products/unconjugated-anti-streptavidin/>