



GOAT ANTI-RABBIT IGG ANTIBODY (H+L), UNCONJUGATED

SKU: AI-1000-1.5



DESCRIPTION

Unconjugated Goat Anti-Rabbit IgG Antibody binds to rabbit primary antibodies for immunohistochemistry and many other applications.

Features:

- Affinity-purified, ultrapure, high affinity antibody
- Thoroughly adsorbed against serum and immunoglobulins from potentially interfering species
- Recognizes both heavy and light chains (H+L)
- Ready for iodination, fluorochrome labeling, or enzyme conjugations
- Can be employed as a capture antibody in enzyme immunoassays or in other assays requiring carrier-free immunoglobulins
- Supplied in solution

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



SPECIFICATIONS

Format	Concentrate
Formulation	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose
Unit Size	1.5 mg
Storage Instructions	2-8 °C; Store frozen for long term storage
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs
Target Species	Rabbit
Concentration	3 mg/ml
Conjugate	Unconjugated
Host Species	Goat

TECHNICAL INFORMATION

The goat anti-rabbit Ig antibodies are prepared by hyperimmunizing animals in a manner that produces high affinity antibodies. These are then purified by an affinity chromatography procedure designed to remove any low affinity antibodies which may be present. Cross-reactivities that are likely to interfere with specific labeling are removed by solid-phase adsorption techniques. The final product is then subjected to rigorous quality control assays including immunodiffusion, solid-phase enzyme immunoassays, gel electrophoresis and solid-phase binding assays. In preparing the labeled antibodies, great care is taken to ensure the maximum degree of labeling with no alteration in the specificity and affinity of the antibody. The labeled antibody then undergoes a further series of quality control assays, including immunohistochemical analysis.

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[Goat Anti-Rabbit IgG Antibody \(H+L\), Unconjugated](https://vectorlabs.com/products/unconjugated-goat-anti-rabbit-igg/)

<https://vectorlabs.com/products/unconjugated-goat-anti-rabbit-igg/>