

www.vectorlabs.com
Email: customerservice@vectorlabs.com
Telephone: (650) 697-3600

MB 680R NHS ESTER

SKU: FP-1671



DESCRIPTION



MB™ 680R is a bright and photostable far-red dye that emits fluorescence at about 685 nm in the borderline spectral region between far-red and near-IR. Although the absorption maximum is at around 680 nm, this dye can be sufficiently excited by the 633 or 635 nm laser. MB™ 680R dye is water soluble and pH-insensitive from pH 4 to pH 10. MB 680R is a rhodamine-based dye, and like rhodamine dyes in general, it is exceptionally photostable. The superior photostability

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

and excellent brightness of MB 680R make the dye an ideal choice for confocal microscopy and other demanding applications.


The NHS Ester reacts specifically and efficiently with a primary amine (e.g., side chain of lysine residues or aminosilane-coated surfaces) at pH 7-9 to form a stable, covalent amide bond. The NHS ester (or succinimidyl ester) is the most popular tool for conjugating dyes to the primary amines of protein or antibody (Lys), amine-modified oligonucleotides, and other amine-containing molecules.

SPECIFICATIONS

Molecular Weight	839.91 (protonated)
Extinction Coefficient	135,000 cm-1M-1
Reactivity	Primary amine
Unit Size	1 mg, 5 mg, 25 mg, 100 mg
Solubility	Water, DMSO, DMF
Storage Instructions	-20°C.
Spectrally Similar Dyes	Alexa Fluor® 680, CF® 680R
Excitation/Emission Maximum	685/709 nm
Shipping Conditions	Ambient temperature
Shipping Instructions	Ambient temperature

ABS/EM SPECTRA


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- Safety Data Sheet
- Datasheet

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.