
www.vectorlabs.com
Email: customerservice@vectorlabs.com
Telephone: (650) 697-3600

NUCLEAR FAST RED COUNTERSTAIN

SKU: H-3403-500



DESCRIPTION

Nuclei stained with Vector® Nuclear Fast Red (H-3403) have a pink appearance. Vector Nuclear Fast Red is an excellent choice in multiple labeling applications when hematoxylin obscures the substrate colors.

Features:

- Good contrast with a variety of substrates
- Rapid staining procedure requiring only one step
- Suitable for use with non-aqueous mountants

SPECIFICATIONS

Format	Ready-to-Use

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

www.vectorlabs.com
Email: customerservice@vectorlabs.com
Telephone: (650) 697-3600

Unit Size 500 ml

Applications Immunohistochemistry / Immunocytochemistry, In situ hybridization

TECHNICAL INFORMATION

Optimize Counterstain Intensity


CITATIONS

Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- User Guide
- Counterstain/Substrate Compatibility Chart
- Safety Data Sheet
- Download CoA
- Datasheet

GALLERY IMAGES





For research use only. Not intended for therapeutic or diagnostic use in animals or humans.